

Cultivation of Vero Cells to High Cell Densities for Human Influenza Vaccine Production

Nandini Prabhakar *, Marcus Wolschek⁺, Micaela Goh Yue Lyn*, Manfred Reiter⁺, Joachim Seipelt⁺, Xiangliang Lin*

*Esco Aster Pte Ltd, Singapore⁺ Nuvonis Technologies, Austria



Introduction

Annual vaccination is an effective method to prevent influenza infection. However, catering to the need of influenza vaccines in developing countries faces many challenges. Timely production and deployment of influenza vaccines is required in non-vaccine producing countries. We focus on the rapid development of influenza vaccine with the goal of generating a perfect bioprocessing solution as an alternative to chicken eggs that can aid to speed up the manufacturing process.

Esco Aster focuses on high-quality, biomanufacturing of vaccines, biologics, and cell-therapy products. To enable our primary vision to help non-vaccine producing countries to be self-sufficient in manufacturing, storing, and distribution of vaccines, we recently partnered with biotech company Nuvonis (Austria) to establish efficient bioprocessing workflows. It will enable the generation of influenza virus using Nuvonis serum-free Vero cell banks. We demonstrate that the Tide Motion manufacturing platform, modularly integrated with Esco Cell Processing Isolator, helps to localize vaccine production, making it more affordable – in terms of CAPEX & OPEX – for developing countries.



Vero Cells

Vero cells are anchorage-dependent cells that are widely used for vaccine production. They are derived from the epithelial kidney cells of African Green monkeys. They have many advantages in terms of high viral infectivity and thus are very effective for primary virus isolation. Vero cells are susceptible to an array of virus such as influenza, rabies, reovirus and Japanese encephalitis virus to name a few. It is desirable to cultivate Vero cells to density compared to those traditionally grown in 2D T-flasks and roller bottles. We culture Vero cells to an optimal density using macrocarriers – BioNOC II in 500ml CelCradle.

Vero Cell Cultivation in Tide Motion Bioreactors

A serum-free Vero cell research cell bank (passage 144-160) is used for all applications. The cell banks – both MCB and WCB – have been fully characterized including tumorigenicity testing at the end of production level (EOP). Serum-free media was supplemented with L-glutamine before use. A serum-supplemented media was also used in cell growth assays for optimization. Glucose in the media was measured using Esco GlucCell™ device and test strips.

Growth kinetics are displayed in Figure 1 (a) and (b) as below. Glucose consumption and Glucose utilization rate are shown in (c) and (d).




Figure 1 Vero cell growth curve in (a) Vero Plus medium and (b) serum-containing medium and the Glucose concentration in the media of the respective cultures. M=media change; HM=half media change; GI=Glucose





Figure 2 Vero cell growth supported by BioNOC II macrocarriers, which provide a large surface area for growth. (L) BioNOC II macrocarriers and (R) Vero cells under 4x magnification. The cells are stained with Fluorescein diacetate (FDA) and Hoechst stain on the extreme right. Left panel represents early culture and right panel represents late culture.

Cultivation of Vero Cells for Influenza Vaccine Production

Cells from five confluent T-175 flasks were harvested by trypsinisation, centrifuged for 5 min at 275 x g, and resuspended in cell growth medium. Cells were mixed with cell growth medium in a total volume of 500 ml and transferred to a CelCradle-500AP (perfusion) using the following Tide Motion parameters. After 1 hour and 45 minutes, media was sampled, and cells were counted using the trypan blue dye exclusion test for cell viability. 85% of the cells attached stably to the carriers. Tide parameters for cell growth are as follows:

Uprate	Uphold	Downrate	Downhold
2.0 mm/s	20 sec	2.0 mm/s	0 sec



Virus Infection Method

A total cell number of 1.7e9 cells was achieved at 167 hours (7 days) after seeding of the Vero cells in serum-free condition. Vero cells were infected with an influenza A strain at multiplicity of infection (MOI) of 0.01. Following viral infection, the medium was replaced seven times. Infectious virus titre in each medium fraction was determined by fluorescent focus assay (FFA). In total 10.45e10 FFU/ml infectious virus was recovered. The schematic is depicted in Figure 3 and the results in Figure 4(a-e).




Figure 3 Vero seed train preparation and Influenza A production process in the CelCradle-500AP; Working Cell Bank (WCB); Multiplicity of Infection (MOI); Fluorescent Focus Assay (FFA).

Virus Yield




Figure 4 (a) CelCradle-500AP: Vero cell growth in serum-free medium (SFM). The red square indicates the time of infection at 161 hours after seeding (b) Glucose concentration in the media (c) Infectious virus titre in log FFU/ml (d) Infectious virus yield in FFU per 500ml fraction (e) High cell density was achieved in 3D culturing compared to 2D culturing

Conclusion

In this proof-of concept study, it has been demonstrated that Vero cells can be grown to high densities culture of 2.9e9 (5e6 cells/ml) using SFM and 3.3e9 (6e6cells/ml) using serum-containing media. This represents an efficient bioprocess workflow of a single-use strategy that can be scaled up for biomanufacturing of affordable vaccines using the Tide Motion platform and Cell Processing Isolator.

Future Studies

These promising results can be further optimized in a cGMP facility, scaling up for the bulk production of influenza vaccines.

Esco Aster Capabilities

Vaccines process developments

Technology transfer of processes worldwide

Acknowledgements

Our acknowledgements and special thanks are due to Nuvonis Technologies GmbH, a biotechnology company based in Vienna Austria, for its fully characterized GMP Vero cell lines and Influenza strains.